

Product Datasheet

PDZ domain-containing protein 11 antibody (FITC) (orb238804)

Description	PDZ domain-containing protein 11 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human PDZ domain-containing protein 11 protein (1-140AA)
Target	PDZD11
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5EBL8
Expiration Date	12 months from date of receipt.