

Product Datasheet

60S ribosomal protein L11 antibody (FITC) (orb238829)

Description	60S ribosomal protein L11 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human 60S ribosomal protein L11 protein (3-178AA)
Target	IgG
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62913
Dilution Range	IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.