

Product Datasheet

Cytokeratin 18 antibody (Biotin) (orb238884)

Description	Cytokeratin 18 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Keratin, type I cytoskeletal 18 protein (102-396AA)
Target	KRT18
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05783
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.