



Product Datasheet

Collagenase 3 antibody (orb238934)


Description	Collagenase 3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Collagenase 3 protein (104-471AA)
Target	MMP13
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P45452
Dilution Range	WB:1:500-1:5000, IHC-P:1:20-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat whole cell...



Western blot analysis of PC3 whole cell ...



Mouse kidney, lung tissue using Collagen...