


Product Datasheet

Gastricsin antibody (orb238946)

Description	Gastricsin antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Gastricsin protein (17-388AA)
Target	PGC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P20142
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse stomach t...