



## Product Datasheet

### Replication protein A 32 kDa subunit antibody (orb238995)


|                            |                                                                          |
|----------------------------|--------------------------------------------------------------------------|
| <b>Description</b>         | Replication protein A 32 kDa subunit antibody                            |
| <b>Species/Host</b>        | Rabbit                                                                   |
| <b>Reactivity</b>          | Human                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                             |
| <b>Tested Applications</b> | ELISA, IF, IHC, IP, WB                                                   |
| <b>Immunogen</b>           | Recombinant Human Replication protein A 32 kDa subunit protein (1-267AA) |
| <b>Target</b>              | RPA2                                                                     |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol           |
| <b>Form/Appearance</b>     | Liquid                                                                   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.            |
| <b>Note</b>                | For research use only                                                    |
| <b>Isotype</b>             | IgG                                                                      |
| <b>Clonality</b>           | Polyclonal                                                               |
| <b>Uniprot ID</b>          | <b>P15927</b>                                                            |
| <b>Dilution Range</b>      | WB:1:1000-1:5000, IHC-P:1:20-1:200, IF:1:50-1:200, IP:1:200-1:2000       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                          |



Western blot analysis of (Lane 1): EC109...



Immunoprecipitation analysis of Hela who...



Immunofluorescence analysis of Hela cell...