

Product Datasheet

Replication protein A 32 kDa subunit antibody (Biotin) (orb238996)

Description Replication protein A 32 kDa subunit antibody (Biotin)

Species/Host Rabbit

Reactivity Human

Conjugation Biotin

Tested Applications ELISA

Immunogen Recombinant Human Replication protein A 32 kDa subunit protein (1-267AA)

Target RPA2

Preservatives Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.


Note For research use only

Isotype IgG

Clonality Polyclonal

Uniprot ID **P15927**

Expiration Date 12 months from date of receipt.



Western blot analysis of (Lane 1): EC109...