

Product Datasheet

X-ray repair cross-complementing protein 6 antibody (orb238999)

Description X-ray repair cross-complementing protein 6 antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC, WB

Immunogen Recombinant Human X-ray repair cross-complementing protein 6 protein (2-609AA)

Target XRCC6

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C.
Avoid repeated freeze.

Note For research use only

Isotype IgG

Clonality Polyclonal

Uniprot ID **P12956**

Dilution Range WB:1:500-1:5000, IHC-P:1:20-1:200,
IF:1:50-1:200

Expiration Date 12 months from date of receipt.

