



Product Datasheet

60S acidic ribosomal protein P1 antibody (orb239012)


Description	60S acidic ribosomal protein P1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human 60S acidic ribosomal protein P1 protein (2-114AA)
Target	RPLP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05386
Dilution Range	WB:1:500-1:2000, IHC-P:1:20-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cells



Immunohistochemical staining of human placenta



Western blot analysis of EC109 whole cell lysate