


Product Datasheet

60S ribosomal protein L17 antibody (orb239016)

Description	60S ribosomal protein L17 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human 60S ribosomal protein L17 protein (2-184AA)
Target	RPL17
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P18621
Expiration Date	12 months from date of receipt.



Western blot analysis of EC109 whole cel...