

Product Datasheet

Zinc finger BED domain-containing protein 1 protein antibody (Biotin) (orb239041)

Description Zinc finger BED domain-containing protein 1
protein antibody (Biotin)

Species/Host Rabbit

Reactivity Human

Conjugation Biotin

Tested Applications ELISA

Immunogen Recombinant Human Zinc finger BED
domain-containing protein 1 protein (1-
210AA)

Target ZBED1

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C. Avoid
repeated freeze.


Note For research use only

Isotype IgG

Clonality Polyclonal

Uniprot ID **096006**

Expiration Date 12 months from date of receipt.



Western
blot
analysis
of (Lane
1):
EC109...