

Product Datasheet

Zinc finger BED domain-containing protein 1 protein antibody (HRP) (orb239043)

Description	Zinc finger BED domain-containing protein 1 protein antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Zinc finger BED domain-containing protein 1 protein (1-210AA)
Target	ZBED1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	096006
Expiration Date	12 months from date of receipt.