


Product Datasheet

39S ribosomal protein L9, mitochondrial antibody (orb239067)

Description	39S ribosomal protein L9, mitochondrial antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human 39S ribosomal protein L9, mitochondrial protein (60-267AA)
Target	MRPL9
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BYD2
Expiration Date	12 months from date of receipt.



Western blot analysis of EC109 (Lane 1), ...