

Product Datasheet

40S ribosomal protein S11 antibody (Biotin) (orb239073)

Description	40S ribosomal protein S11 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human 40S ribosomal protein S11 protein (2-158AA)
Target	RPS11
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62280
Expiration Date	12 months from date of receipt.