


Product Datasheet

40S ribosomal protein S9 antibody (orb239091)

Description	40S ribosomal protein S9 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human 40S ribosomal protein S9 protein (2-184AA)
Target	RPS9
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P46781
Expiration Date	12 months from date of receipt.



Western blot analysis of EC109 (Lane 1), ...