



Product Datasheet

6.8 kDa mitochondrial proteolipid antibody (orb239095)

Description	6.8 kDa mitochondrial proteolipid antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human 6.8 kDa mitochondrial proteolipid protein (1-58AA)
Target	MP68
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P56378
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of MCF-7
cel...