

Product Datasheet

60 kDa SS-A/Ro ribonucleoprotein antibody (orb239099)

Description 60 kDa SS-A/Ro ribonucleoprotein antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC, WB

Immunogen Recombinant Human 60 kDa SS-A/Ro ribonucleoprotein protein (3-535AA)

Target TROVE2

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Note For research use only


Isotype IgG

Clonality Polyclonal


Uniprot ID **P10155**

Dilution Range WB:1:500-1:5000, IHC-P:1:20-1:200, IF:1:50-1:200


Expiration Date 12 months from date of receipt.



Western blot analysis of Hela (Lane 1), ...



Immunohistochemical staining of human pr...



Western blot analysis of Hela whole cell...