



Product Datasheet

Bis(5'-adenosyl)-triphosphatase antibody (orb239147)

Description	Bis(5'-adenosyl)-triphosphatase antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Bis(5'-adenosyl)-triphosphatase protein (2-147AA)
Target	FHIT
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49789
Dilution Range	IHC-P:1:500-1:1000, IF:1:50-1:500
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HepG2 cells.



Immunohistochemistry of paraffin-embedded tissue.