

Product Datasheet

DNA (cytosine-5) methyltransferase3A antibody (HRP) (orb239166)

Description	DNA (cytosine-5) methyltransferase3A antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human DNA (cytosine-5)-methyltransferase 3A protein (680-902AA)
Target	DNMT3A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y6K1
Expiration Date	12 months from date of receipt.