

Product Datasheet

DNA repair protein complementing XP-C cells antibody (FITC) (orb239169)

Description	DNA repair protein complementing XP-C cells antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human DNA repair protein complementing XP-C cells protein (496-734AA)
Target	XPC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q01831
Expiration Date	12 months from date of receipt.