

Product Datasheet

Histone H2A.J antibody (Biotin) (orb239224)

Description	Histone H2A.J antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Histone H2A.J protein (2-129AA)
Target	H2AFJ
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BTM1
Expiration Date	12 months from date of receipt.