

Product Datasheet

papillomavirus type 16 Protein E6 antibody (orb239255)

Description	papillomavirus type 16 Protein E6 antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human papillomavirus type 16 Protein E6 protein (1-158AA)
Target	E6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03126
Expiration Date	12 months from date of receipt.