


Product Datasheet

Prefoldin subunit 5 antibody (orb239265)

Description	Prefoldin subunit 5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Prefoldin subunit 5 protein (2-154AA)
Target	PFDN5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99471
Expiration Date	12 months from date of receipt.



Western blot analysis of EC109(Lane 1), ...