



## Product Datasheet

### Protein arginine N-methyltransferase 1 antibody (orb239276)


|                            |                                                                            |
|----------------------------|----------------------------------------------------------------------------|
| <b>Description</b>         | Protein arginine N-methyltransferase 1 antibody                            |
| <b>Species/Host</b>        | Rabbit                                                                     |
| <b>Reactivity</b>          | Human, Mouse                                                               |
| <b>Conjugation</b>         | Unconjugated                                                               |
| <b>Tested Applications</b> | ELISA, IF, IHC, WB                                                         |
| <b>Immunogen</b>           | Recombinant Human Protein arginine N-methyltransferase 1 protein (1-248AA) |
| <b>Target</b>              | PRMT1                                                                      |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol             |
| <b>Form/Appearance</b>     | Liquid                                                                     |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.              |
| <b>Note</b>                | For research use only                                                      |
| <b>Isotype</b>             | IgG                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                 |
| <b>Uniprot ID</b>          | <b>Q99873</b>                                                              |
| <b>Dilution Range</b>      | WB:1:500-1:5000, IHC-P:1:20-1:200, IF:1:50-1:200                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                            |



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of Mouse kidney t...



Immunofluorescence staining of A549 cell...