



Product Datasheet

Signal recognition particle 19 kDa protein antibody (orb239312)


Description	Signal recognition particle 19 kDa protein antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Signal recognition particle 19 kDa protein (2-144AA)
Target	SRP19
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09132
Dilution Range	WB:1:1000-1:5000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1):jurkat...



Immunohistochemical staining of human pa...



Immunohistochemical staining of human li...