


Product Datasheet

SOSS complex subunit B2 antibody (orb239316)

Description	SOSS complex subunit B2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human SOSS complex subunit B2 protein (1-204AA)
Target	NABP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96AH0
Expiration Date	12 months from date of receipt.



Western blot analysis of jurkat(Lane 1),...