

Product Datasheet

Ubiquitin-conjugating enzyme E2 variant 2 antibody (orb239336)

Description Ubiquitin-conjugating enzyme E2 variant 2 antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IHC

Immunogen Recombinant Human Ubiquitin-conjugating enzyme E2 variant 2 protein (2-145AA)

Target UBE2V2

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C.
Avoid repeated freeze.

Note For research use only


Isotype IgG

Clonality Polyclonal

Uniprot ID **Q15819**

Dilution Range IHC-P:1:20-1:200

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of human br...