


Product Datasheet

Vesicular, overexpressed in cancer, prosurvival protein 1 antibody (orb239340)

Description	Vesicular, overexpressed in cancer, prosurvival protein 1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Vesicular, overexpressed in cancer, prosurvival protein 1 protein (23-172AA)
Target	VOPP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96AW1
Expiration Date	12 months from date of receipt.



Western blot analysis of EC109(Lane 1), ...