

Product Datasheet

Peptidyl-prolyl cis-trans isomerase A antibody (orb239388)

Description Peptidyl-prolyl cis-trans isomerase A antibody

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications ELISA, WB

Immunogen Recombinant Human Peptidyl-prolyl cis-trans isomerase A protein (1-165AA)

Target PPIA

Preservatives Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.


Note For research use only

Isotype IgG


Clonality Polyclonal

Uniprot ID **P62937**

Expiration Date 12 months from date of receipt.



Western blot analysis of HepG2 cell, Mo...



Western blot analysis of 293T whole cel...