

Product Datasheet

Advanced glycosylation end product-specific R antibody (Biotin) (orb239397)

Description	Advanced glycosylation end product-specific R antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Advanced glycosylation end product-specific receptor protein (23-342AA)
Target	AGER
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15109
Expiration Date	12 months from date of receipt.