

Product Datasheet

C-C motif chemokine 14 antibody (FITC) (orb239423)

Description	C-C motif chemokine 14 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human C-C motif chemokine 14 protein (20-93AA)
Target	CCL14
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16627
Expiration Date	12 months from date of receipt.