



Product Datasheet

Complement alpha 2 Macroglobulin antibody (orb239453)


Description	Complement alpha 2 Macroglobulin antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Complement C3 protein (26-225AA)
Target	C3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01024
Dilution Range	IHC-P:1:50-1:500, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pr...



Immunohistochemical staining of human ce...



Immunohistochemical staining of human br...