


Product Datasheet

C-X-C motif chemokine 13 antibody (orb239465)

Description	C-X-C motif chemokine 13 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human C-X-C motif chemokine 13 protein (23-95AA)
Target	CXCL13
Preservatives	0.1mol/L NaCl, 15mmol/L NaN ₃ , PH7.2
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	043927
Dilution Range	IHC: 1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human li...