

Product Datasheet

C-X-C motif chemokine 5 antibody (Biotin) (orb239470)

Description	C-X-C motif chemokine 5 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human C-X-C motif chemokine 5 protein (37-110AA)
Target	CXCL5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P42830
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T whole cell...