

Product Datasheet

Granulocyte-macrophage colony-stimulating factor antibody (Biotin) (orb239490)

Description	Granulocyte-macrophage colony-stimulating factor antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Granulocyte-macrophage colony-stimulating factor protein (18-144AA)
Target	CSF2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04141
Expiration Date	12 months from date of receipt.