



Product Datasheet

Granzyme B antibody (orb239493)


Description	Granzyme B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Granzyme B protein (21-247AA)
Target	GZMB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10144
Dilution Range	IHC-P:1:500-1:1000, IF:1:200-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse stomach t...



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of human sp...