



Product Datasheet

Hepatitis A virus cellular R 2 antibody (orb239497)


Description	Hepatitis A virus cellular R 2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Hepatitis A virus cellular receptor 2 protein (22-202AA)
Target	HAVCR2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TDQ0
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of using Hep...



Immunofluorescence analysis of MCF-7 cel...



Immunohistochemical staining of Mouse ki...