



Product Datasheet

IgG R FcRn large subunit p51 antibody (orb239550)


Description	IgG R FcRn large subunit p51 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human IgG receptor FcRn large subunit p51 protein (24-297AA)
Target	FCGRT
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P55899
Dilution Range	WB:1:500-1:2000, IHC-P:1:20-1:1000, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis
of HepG2 whole cel...



Immunohistochemical
staining of human
pl...



Immunohistochemical
staining of human
Pa...