

## Product Datasheet

### IgG R FcRn large subunit p51 antibody (FITC) (orb239552)

|                        |                                                                          |
|------------------------|--------------------------------------------------------------------------|
| <b>Description</b>     | IgG R FcRn large subunit p51 antibody (FITC)                             |
| <b>Species/Host</b>    | Rabbit                                                                   |
| <b>Reactivity</b>      | Human                                                                    |
| <b>Conjugation</b>     | FITC                                                                     |
| <b>Immunogen</b>       | Recombinant Human IgG receptor FcRn large subunit p51 protein (24-297AA) |
| <b>Target</b>          | FCGRT                                                                    |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol           |
| <b>Form/Appearance</b> | Liquid                                                                   |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C.<br>Avoid repeated freeze.         |
| <b>Note</b>            | For research use only                                                    |
| <b>Isotype</b>         | IgG                                                                      |
| <b>Clonality</b>       | Polyclonal                                                               |
| <b>Uniprot ID</b>      | <b>P55899</b>                                                            |
| <b>Expiration Date</b> | 12 months from date of receipt.                                          |