

Product Datasheet

Killer cell lectin-like R subfamily G member 1 antibody (orb239630)

Description	Killer cell lectin-like R subfamily G member 1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Killer cell lectin-like receptor subfamily G member 1 protein (60-195AA)
Target	KLRG1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96E93
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.