

Product Datasheet

Non-secretory ribonuclease antibody (HRP) (orb239649)

Description	Non-secretory ribonuclease antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Non-secretory ribonuclease protein (28-161AA)
Target	RNASE2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10153
Expiration Date	12 months from date of receipt.