



Product Datasheet

Protein-glutamine gamma-glutamyltransferase 2 antibody (orb239690)


Description	Protein-glutamine gamma-glutamyltransferase 2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Protein-glutamine gamma-glutamyltransferase 2 protein (4-544AA)
Target	TGM2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P21980
Dilution Range	WB:1:1000-1:5000, IHC-P:1:200-1:500, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 whole cel...



Immunofluorescence analysis of Hela cell...



Immunohistochemical staining of human li...