

Product Datasheet

Src kinase-associated phosphoprotein 1 antibody (orb239702)

Description	Src kinase-associated phosphoprotein 1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Src kinase-associated phosphoprotein 1 protein (1-358AA)
Target	SKAP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86WV1
Expiration Date	12 months from date of receipt.