



Product Datasheet

Transcription initiation factor IIB antibody (orb239738)


Description	Transcription initiation factor IIB antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Transcription initiation factor IIB protein (1-316AA)
Target	GTF2B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q00403
Dilution Range	WB:1:500-1:5000, IHC-P:1:20-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 (Lane 1), ...



Immunohistochemical staining of human kidney tissue.



Immunofluorescence analysis of PC-3 cells.