



Product Datasheet

CD71 antibody (orb239742)


Description	CD71 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Transferrin receptor protein 1 protein (94-754AA)
Target	TFRC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02786
Dilution Range	IHC-P: 1:20-1:200, IF/ICC: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human CO...



Immunofluorescence analysis of HepG2 cel...



Immunohistochemical staining of human CO...