


Product Datasheet

Tumor necrosis factor R superfamily member 11A antibody (Biotin) (orb239752)

Description	Tumor necrosis factor R superfamily member 11A antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 11A protein (28-202AA)
Target	TNFRSF11A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y6Q6
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1): EC109...