

Product Datasheet

Alcohol dehydrogenase class 4 mu/sigma chain antibody (orb239832)

Description Alcohol dehydrogenase class 4 mu/sigma chain antibody

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC, WB

Immunogen Recombinant Human Alcohol dehydrogenase class 4 mu/sigma chain protein (1-386AA)

Target Ribosome-inactivating protein karasurin-C

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C.
Avoid repeated freeze.

Note For research use only


Isotype IgG

Clonality Polyclonal


Uniprot ID **P40394**

Dilution Range WB:1:1000-5000, IHC-P:1:20-200


Expiration Date 12 months from date of receipt.



Western blot analysis of Mouse liver tis...



Immunohistochemical staining of human lu...



Immunohistochemical staining of human ga...