

Product Datasheet

Alcohol dehydrogenase class 4 mu/sigma chain antibody (Biotin) (orb239833)

Description	Alcohol dehydrogenase class 4 mu/sigma chain antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Alcohol dehydrogenase class 4 mu/sigma chain protein (1-386AA)
Target	Botulinum neurotoxin type E
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P40394
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.