

Product Datasheet

Cholesterol side-chain cleavage enzyme, mitochondrial antibody (orb239864)

Description	Cholesterol side-chain cleavage enzyme, mitochondrial antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Cholesterol side-chain cleavage enzyme, mitochondrial protein (40-521AA)
Target	CYP11A1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05108
Dilution Range	WB:1:1000-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of rat adrenal gla...