

Product Datasheet

Citrate synthase, mitochondrial antibody (orb239868)

Description Citrate synthase, mitochondrial antibody

Species/Host Rabbit

Reactivity Human, Zebrafish

Conjugation Unconjugated

Tested Applications ELISA, IHC, WB

Immunogen Recombinant Human Citrate synthase, mitochondrial protein (28-464AA)

Target CS

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Note For research use only


Isotype IgG

Clonality Polyclonal


Uniprot ID **075390**

Dilution Range WB:1:500-5000, IHC-P:1:20-200


Expiration Date 12 months from date of receipt.



Western blot analysis of HeLa (Lane 1), H...



mitochondrial antibody



Immunohistochemical staining of human pr...